## LOCATORS



#### LT-490 Locator

Receives at 433.92Mhz from wireless transmitters and retransmits at 922 MHz to either a repeater (RP-990) and or master receiver (MR 500). Adjustable receive range for more granular locate capability.

**MR-500** 



## RcareInc.com

#### **RP-990** Repeater

The repeater is designed to increase the reach of the locate system by increasing the distance between control and locators. It receives 922MHz signals from both LT-490 locators and RP-990 repeaters.

#### **Heated and Cooled Enclosure**

#### WPE HC-1





14" x 12" x 6" heated/cooled NEMA4 locking latch enclosure. Includes: heater, fan, thermostats, strain relief, 2 bulkhead holes, screws and 2 power plugs. Durable polycarbonate for use with locators or repeaters in harsh outdoor environments such as parking lots, walk paths or golf courses.



Response Care is your provider of highly reliable, easy to use equipment and software for the assisted and independent living industry. Our components are designed to integrate into a variety of communication systems in order to facilitate efficient and verifiable response to medical emergencies.

R-care systems and components are only available from trained, authorized distributors. For an authorized distributor in your area, please call us toll free at 1-888-724-1911.

## IN-ROOM COMMUNICATORS, VOICE TO VOICE



This 3 button unit call, cancel and programmable concierge is a low profile wall or desk mount communicator. Two phone number capability, primary and back-up. For use with voice to voice, on-site monitoring. 9 VDC power with integrated battery back up.



#### **Wall Mount Voice Dialer Communicator** RC-1800

Integrated voice record capability with up to six phone number dial out capability. LCD display with clock, program prompts and status indicators. After an emergency message (up to 20 seconds) care giver can carry on high quality, full duplex conversation with protected resident.

# plastic case.





#### Wireless Personal Emergency For Active Senior Comm

474 Ridge Road, Webster, NY 14580 Tel: 585-671-4144 Fax: 585-671-1183 Info@RcareInc.com

Toll Free: 1-888-724-1911



### RcareInc.com

## **PRODUCT CATALOG**



- In-Room Communicators, Receivers/Dialers Transmitters
- Personal Transmitters Locators Repeaters
- Heated & Cooled Enclosures Software

#### **Apartment Communicator RC-2200** with Concierge Button

### **RC-1100**

Integrated voice record dialer for desk or table with optional wall mount bracket. Many of the features of the RC-1800 but with larger, easy to see, red activation button. Features a long-lasting, high-impact

## **Reare** Wireless Personal Emergency

MS-6

**WM-8** 

## **Door Monitoring for Residential Safety**



#### **Remote Keypad RK-77**

**Automated Check-in** 

**Motion Sensor** 

Wall Mount Push

**Button Transmitter** 

16 button backlit keypad with power conservation and tamper. Used in conjunction with door and window contact modules. Allows authorized use of monitored door. Power: 3 AAA batteries (included).

Infrared motion sensor for use with R-Care

Check-in resident check-in system. Surface

and corner mounting with tamper for

automatic power saver mechanism.

removal. Low Battery reporting, 6 hour

Detection range 12 meters, insect and

draft proof. LED fault and test indicator. Power: 2 AA batteries (included).

Palm-sized design with large red button

common areas for alerts or as a remote

and call assurance LED. For use in

## **Door and Window Contact**

Compact design for easy mounting on door frame. Includes tamper protection, low battery reporting, built in auxiliary input terminal. Built in Reed Switch and powerful magnet gives 30mm of space between door and switch. External contact can be used simultaneously. LED fault and test indicator. Power: 2 AA batteries (included).

WD-3

#### **Bathroom Pull Cord** BP-7R

Bathroom emergency pull cord mounts to a plastic single gang standard electrical box (not included) or ring. Locking call switch connected to a 6 ft. nylon cord. Re-transmit activation until reset. Reports reset and low battery condition. The only wireless pull cord with red LED call assurance indicator. Power: 2 AAA batteries (included).

#### **Universal Transmitter**

#### **UT-3 & UT-3R**

NC-66

Connects closure transmits alarm. LED indicator for activation. Periodic transmission of alarm for power conservation. PCT/R provides automatic low battery warning and transmits restore signal.

## **Personal Transmitters**

reset button.

### Wearable transmitter/Convertible WTC

pendant or convertible to a wristband. Permanently sealed, battery operated and water proof up to 20 ft. Low battery warning.

## **Nurse Call Adapter**

Connects to existing hardwired nurse call systems to provide wireless activation. Up to twenty (20) wearable transmitters can be programmed. Momentary and maintain contact. Battery backup.

## For Active Senior Communities ReareInc.com


## SYSTEMS AND SOFTWARE



#### **R-care Plus Software RC-PLUS VOICE TO VOICE**



R-care PLUS, mission critical computer with intuitive user interface, the ability to deliver resident "alerts" and "voice to voice" to care givers facility wide with telephone and wireless telephone systems. Enhance the security of your facility by monitoring doors, windows, fixed location wireless or wireless pull cords, in "risk areas", as well as environment monitoring. Two-way voice in-room resident communicators utilize the same wearable transmitters that activate R-care Locator. Extensive logging and reporting of all incidents and can even record emergency conversations between staff and residents.









## **NEW, R-cube Monitoring System** An Exponential Leap In Resident-Alert Communications

Real time, multipoint information retrieval is finally here! Receive resident alerts from any device, anytime, anywhere in your facility- It's communications to the third power. Browser based with hardware and software in one, self-contained unit, R-cube provides the ultimate in information retrieval capability and flexibility.